

- 👂 Has your hearing aid stopped working ?
- 👂 Have you checked your battery?
- 👂 Have you checked the tube for condensation?

When the weather gets colder, droplets of water can collect in the plastic tube and 'prevent' the hearing aid from working effectively.

If there is condensation in the tube, the hearing aid will appear as though it is not working.

Pull the tube from the elbow of the hearing aid as though you were going to wash the mould and then either blow through the tube or use a 'puffer' to blow air through.



PUFFER

If problems persist, then contact Hearing Services to arrange an appointment to have the hearing aid checked.

'Puffers' available to purchase from Reception at Hearing Services. Please ask reception or audiologist for latest charges.

Hearing Services Department
West Park Rehabilitation Hospital
Park Road West
Wolverhampton
WV1 4PW

Telephone: 01902 444055

If you require this leaflet in a different format i.e. larger print, please contact:



Telephone: 01902 444055

Available from 8am-17:00pm



Fax: 01902 444056



e mail: rwh-tr.hearing-services@nhs.net

If you require this leaflet in a different format i.e. larger print etc please contact us above.

Approved Date: April 2013
RATIFICATION DATE: Dec 2013
Review Date: April 2018

Reference Number: PIL 13



CONDENSATION IN MY TUBE!

WHAT IS THIS?

SELF HELP LEAFLET